

# BD-DP Disposable video laryngoscope

Anti-fog medical-grade polycarbonat



Neonate Handle

Adult Handle

Item Number:BD-DP

Handle Order No.	Description
BD-DP-H1	Neonate Handle for Miller 0 blades
BD-DP-H2	Adult Handle for MAC1, MAC2, MAC3,MAC4 blades
BD-DP-H3	Adult Handle for D blades

## Blade parameter size

Blades Order No.	DP-Mil0	DP-MAC1	DP-MAC2	DP-MAC3	DP-MAC4	DP-XD
Size	Mil0	MAC1	MAC2	MAC3	MAC4	D Blade
Application	Neonate	Pediatric	Small Adult	Adult	Large Adult	Difficult airway
Length(cm)	8.31	10.15	12.11	13.69	16.08	17.5

## Monitor Specification

	Description	Technical Index
Parts	Display	3inch LCD
	Resolution	640*480 (RGB)
	Viewing angle	≥60°
	Power	<2W
Monitor	Display front and back rotation Angle	0°~130°
	Weight	225g
Camera	Illuminance	≥150 LUX
	Resolution	1280*720 px
Image/Video function	Image/ Video function	Yes
	Output	USB output, Easy to storage and establish image/video
	Memory type	Micro SD card 32GB
Batterty	Battery Type	Rechargeable lithium battery
	Capacity	1350mAh
	Cycle life of battery	>500 times
	Battery working time	>240min
	Charging time	<2(hrs)