

Flexible Videoscope

- Versatile
- Specialized:
bedside diagnosis & treatment





READY FOR BEDSIDE

- lightweight ($\leq 500\text{g}$)
- self-contained

READY FOR EMS

- integrated battery

DURABLE

- heat-resistant, corrosion-resistant, wear-resistant, and non-conductive material
- crystal sapphire screen and camera lens

Insertion tube outer diameter (mm)	Working channel diameter (mm)	Recommended for
Ø 2.6	--	• inspection
Ø 2.8	--	• pediatric airway management
Ø 3.8	--	• check the placement of DLT
Ø 3.1	Ø 1.2	• check the placement of endobronchial blocker
Ø 3.9	Ø 1.5	• endotracheal drug administration & inspection
Ø 4.5	Ø 1.5	• tracheal intubation
Ø 5.2	Ø 2.2	• check the placement of DLT
Ø 5.2	Ø 2.6	• suction sputum
Ø 5.5	Ø 2.4	• BAL procedures
		• remove foreign objects
		• sample lung issue

SPECIFICATIONS

Display	3" TFT LCD	Camera	CMOS
Tilt and Rotation	0°-180° tilt, 0°-180° rotate	Memory	32G
Distal Dip Angle	upward 160°+, downward 130°+	Battery Life	≥ 180 mins on a full charge
Angle of View	90°	Charging Time	< 4 hrs
Depth of field	3-100mm		