



Video Laryngoscope with Single-use Sleeve

UED-A

- **Low cost**
- **Low cross-contamination risk**



UED-M3-S adds on:



TOUCH SCREEN

image mirroring: flips upside-down whenever needed
auto shutdown: powers off after 10 mins of inactivity



WIFI READY

real-time screen sharing
fast access to photos and videos

2 monitor options (UED-M3, UED-M3-S)



UED-A0
newborn, 2.7-3.8kg



UED-A1
0-2 yrs, 3.4-13.2kg



UED-MA1
0-2 yrs, 3.4-13.2kg



UED-A2
2-6 yrs, 13.2-23.6kg



UED-A3
6-adult, above 23.6kg



UED-A4
adult, obesity



Compatible with UE M8 8" monitor

REUSABLE HANDLE

- easy to control

GRAB AND GO

- quick set up
- go-to device for prehospital (EMS) & hospital

SIGNATURE SLEEVE ANGLE

- clearer view of the glottis
- easier tracheal intubation

SPECIFICATIONS

MONITOR		SLEEVES	
Display	3" TFT LCD, 720x480 px	Angle of View	60°
Tilt and Rotation	0°-110° tilt, 0°-270° rotate	Light Source	high intensity LED
Memory	32G, take > 400k photos, record ≥16 hrs of video	Charging Time	< 3 hrs
Battery Life	> 3.5 hrs on a full charge		